

Monoclonal Anti LRH-1 (Liver receptor homolog-1)

Ref: AO-01-LRH-2E1

Volume: 10, 25, 50 and 100 µl

Clone: 1LRH-2E1

Immunogen: Synthetic peptide coupled to Ovalbumin (TMPPHSSLHGYQPYGHF-C) corresponding to aa 267 to 283 of mouse LRH-1

Murine Isotype: IgG2, λ

Specificity : mouse LRH-1

Applications :

- ELISA
- Western Blot
- Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

