

Monoclonal anti methyl K115 Histone H3

Ref: AO-01-H3-3F11

Volume: 10, 25, 50 and 100 µl

Clone: 11HH3-3F11

Immunogen: synthetic peptide coupled to ovalbumin (C-AIHA[me K]RVTIM) corresponding to aa 111-120 of human Histone H3. This sequence is conserved in a wide range of species.

Murine Isotype: IgG2a, K

Specificity: methylated K115 Histone H3 (17 kDa).

Applications :

- ELISA
- Western blot
- Immunofluorescence

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years). Avoid repeat freezing and thawing cycles.

Does not contain any preservative