

Anti- β -Tubulin mAb

Ref: AO-01-TUB-2A2

Volume: 10, 25, 50 and 100 μ l

Clone: 1TUB-2A2

Immunogen: synthetic peptide corresponding to amino acids 412-430 of human β -tubulin (CESNMNDLVSEYQQYQDATA).

Murine Isotype: IgG1, κ

Specificity: human, mouse, rat, ...

Applications :

- ELISA
- Western blot
- Immunoprecipitation
- Immunofluorescence

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years) Avoid repeat freezing and thawing cycles.

Does not contain any preservative

