

## Anti-Actin mAb

**Ref:** AO-01-ACT-2D7

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 1ACT-2D7

**Immunogen:** synthetic peptide corresponding to amino acids 364-375 of human actin (CYDEAGPSIVHRK).

**Murine Isotype:** IgG2a, κ

**Specificity:** human, mouse, rat, ...

### **Applications :**

- ELISA
- Western blot
- Immunocytochemistry
- Immunofluorescence

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years) Avoid repeat freezing and thawing cycles.

***Does not contain any preservative***