

Anti-CETP (CHOLESTERYL ESTER TRANSFER PROTEIN) mAb

Ref: AO-01- CETP-2C5

Volume: 10, 25, 50 and 100 µl

Clone: 1CET-2C5

Immunogen: synthetic peptide coupled to ovalbumin: aa 143-160 of human CETP: NTQLTCDSGRVRTDAPDC.

Murine Isotype: IgG1, κ

Specificity: human CETP.

Applications :

- ELISA
- Western Blot
- Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years).

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

