

## Anti-CETP (CHOLESTERYL ESTER TRANSFER PROTEIN) mAb

**Ref:** AO-01-CETP-2H6

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 2CET-2H6

**Immunogen:** synthetic peptide coupled to ovalbumin: aa 482-493 of human CETP: EHLLVDFLQSLs.

**Murine Isotype:** IgG1, κ

**Specificity:** human CETP.

**Applications :**

- ELISA
- Western Blot
- Immunocytochemistry

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years).

Avoid repeat freezing and thawing cycles.

*Does not contain any preservative*

