

Monoclonal Anti-EMI-2 (Endogenous meiotic inhibitor 2/F-box only protein 43)

Ref: AO-01-EMI-2A10

Volume: 10, 25, 50 and 100 µl

Clone: 2EMI-2A10

Immunogen: synthetic peptide: CRRGSAKARGSKDVL corresponding to aa 613 to 627 of mouse EMI-2.

Murine Isotype: IgA, κ

Specificity : mouse EMI-2.

Applications :

- ELISA
- Western Blot
- Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

