

## Anti-FLAG mAb

**Ref:** AO-01- 2EL-1B11

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 2EL-1B11

**Immunogen:** Synthetic peptide coupled to ovalbumin MDYKDDDDK(C)

**Murine Isotype:** IgG1, κ

**Specificity:** human, mouse, rat, ...

### Applications :

- ELISA
- Western blot
- Immunoprecipitation
- Immunofluorescence

**Recommended dilutions:** 1/1000-1/5000

**Storage:** Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

*Does not contain any preservative*

