

## **Anti-FAP (Fibroblast Activation Protein) mAb**

**Ref:** AO-01-FAP-1E1

**Volume:** 10, 25, 50 and 100 µl

**Clone:** MYC-1A1

**Immunogen:** Synthetic peptide coupled to Ovalbumin (NNIYLKQRPGDPPFQIT-C) corresponding to aa 168 to 184 of mouse FAP.

**Murine Isotype:** IgG1, κ

**Specificity:** mouse and human FAP

**Applications :**

- ELISA
- Western Blot
- Immunofluorescence
- 

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

*Does not contain any preservative*

