

Anti-GFP mAb

Ref: AO-01- GFP-2A3

Volume: 10, 25, 50 and 100 μ l

Clone: 1GFP-2A3

Immunogen: eGFP protein

Murine Isotype: IgG1, κ

Specificity: GFP, GFP-fusion proteins

Applications :

- ELISA
- Western blot
- Immunoprecipitation
- Immunofluorescence

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

