

Monoclonal Anti-GLUTATHION-S-TRANSFERASE Anti-GST

Ref: AO-01-GST-1D10

Volume: 10, 25, 50 and 100 µl

Clone: 15TF2-1D10

Immunogen: GST fusion protein (pGEX Vector)

Murine Isotype: IgG1, κ

Specificity: GST and GST fusion protein pGEX 2T and 3X

Applications :

- Immunoblotting
- ELISA
- Immunoprecipitation
- Immunocytochemistry : on transfected cells
- ChIP
-

Recommended dilutions: 1/1000-1/5000

Storage: Store at -20°C (up to 3 years). Avoid repeat freezing and thawing cycles.

Does not contain any preservative

Reference:

Yue Zhao, Guillaume Lang, Saya Ito, Jacques Bonnet, Eric Metzger, Shun Sawatsubashi, Eriko Suzuki, Xavier Le Guezennec, Hendrik G. Stunnenberg, Aleksey Krasnov, Sofia G. Georgieva, Roland Schüle, Ken-Ichi Takeyama, Shigeaki Kato, László Tora, Didier Devys. A TFTC/STAGA Module Mediates Histone H2A and H2B Deubiquitination, Coactivates Nuclear Receptors, and Counteracts Heterochromatin Silencing. Mol Cell 2008. 29: 92–101