

Anti-Gal 4 (DNA Binding Domain) mAb

Ref: AO-01-GV-3

Volume: 10, 25, 50 and 100 µl

Clone: 2GV-3

Immunogen: Gal 4 DBD fusion protein from aa 1 to aa 147

Murine Isotype: IgG1, κ

Specificity: Gal 4 DBD containing proteins, epitope located in the second half of the DBD from aa 74 to 147

Applications :

- ELISA
- Immunoblotting
- Immunocytochemistry on transfected cells
- Immunoprecipitation
- Supershift experiments
- Transcription assays : see reference

Recommended dilutions: 1/1000-1/5000

Storage: Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

Reference: White & al, EMBO J.(1992) vol.11, n°6 : 2229-2240



