

Monoclonal Anti-HSP60 (Heat Shock Protein 60, Chaperonin, GroEL Homolog)

Ref: AO-01-HSP-2H7

Volume: 10, 25, 50 and 100 µl

Clone: 4MTE-2H7

Immunogen: Natural immunization (unknown immunogen).

Murine Isotype: IgG1, κ

Specificity : Mouse and human HSP 60. Epitope localized in the C-terminal part (190 aa) of human HSP 60

Applications :

- ELISA
- Western Blot
- Immunoprecipitation
- Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

