

## **Anti-His-Tag mAb**

**Ref:** AO-01-HIS-1G4

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 2AMY-1G4

**Immunogen:** 6 Histidine fusion protein.

**Murine Isotype:** IgG1, κ

**Specificity:** 6-Histidine fusion protein (6 His at Cter or Nter)

**Applications :**

- ELISA
- Western blot
- Immunocytochemistry

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

*Does not contain any preservative*

