

Monoclonal Anti-Histone H4

Ref: AO-01- H4-4G8

Volume: 10, 25, 50 and 100 μ l

Format: Crude ascites fluid

Clone: 3HH4-3H10

Immunogen: synthetic peptide: MSGRGKGGKGLGKGGAKRHRKVLK-(C) corresponding to aa 1 to 24 of human histone H4.

This sequence is conserved in a wide range of species

Murine Isotype: IgG1, κ

Specificity: histone H4.

Applications :

- ELISA
- Western blot
- Immunoprecipitation

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years) Avoid repeat freezing and thawing cycles.

Does not contain any preservative