

Anti-Human Huntingtin mAb

Ref: AO-01-HU-2E8

Volume: 10, 25, 50 and 100 µl

Clone: 2HU-2E8

Immunogen: Huntingtin fragment from aa 2146 to aa 2541 as a fusion protein

Murine Isotype: IgG1, κ

Specificity: Huntingtin from human and monkey, reacts weakly or not with mouse protein. Other species unknown. Cross-reaction on other proteins: not detectable by Western blot.

Applications :

- ELISA
- Western blot
- Immunohistochemistry on frozen sections and microwave oven treated paraffin sections (human tissues).
- Immunocytochemistry on transfected cells
- Immunoprecipitation : not tested

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

References :

Trottier & al ; Nature Genetics (1995) 10: 104-110

Trottier & al ; Nature, (1995) 378:403-406



