

## Anti-Human Huntingtin mAb

**Ref:** AO-01-HU-4E6

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 3HU-4E6

**Immunogen:** Huntingtin fragment from aa 1247 to aa 1646 as a fusion protein

**Murine Isotype:** IgG2b, κ

**Specificity:** Huntingtin from human and mouse. Other species unknow. Cross-reaction with other proteins : not detectable by western blot.

### Applications :

- ELISA
- Western blot: not suitable for the natural Huntingtin.
- Immunohistochemistry on frozen sections and microwave oven treated paraffin sections (human tissues).
- Immunocytochemistry on transfected cells
- Immunoprecipitation

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

*Does not contain any preservative*

### References :

Trottier & al ; Nature Genetics (1995) 10: 104-110

Trottier & al ; Nature, (1995) 378:403



