

Anti-Interleukin IL-17A mAb

Ref: AO-01-IL-2B10

Volume: 10, 25, 50 and 100 µl

Clone: 3IL-2B10

Format : Protein G purified

Concentration: 1mg/ml

Buffer: 0.1 M Tris-glycine, pH 7.4.

Immunogen: Ovalbumin conjugated synthetic peptide: LHRNEDPDRYPSVIWEAQC corresponding to aa 79-97 of mouse IL-17A. This sequence has 89% identity with human IL 17A.

Murine Isotype: IgG1, κ

Applications :

- ELISA
- Western Blot
- Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Storage: Store at 4°C

