

Monoclonal Anti-Keratin 5

Ref: AO-01-KER-1A8

Volume: 10, 25, 50 and 100 µl

Clone: 2KR-1A8

Immunogen: Synthetic peptide (FVSTTSSSRKSFKS) corresponding to aa 577-590 of human keratin 5.

Murine Isotype: IgG1, κ

Specificity : mouse and human keratin 5.

Applications :

- ELISA
- Western Blot
- Immunocytochemistry
- Immunohistochemistry

Recommended dilutions: 1/500-1/5000

Storage : store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

