

## Anti-Keratin 5 mAb

**Ref:** AO-01- KER-1C11

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 1KR-1C11

**Immunogen:** Synthetic peptide (MSRQSSVSFRSGGS) corresponding to aa 1-14 of mouse and human keratin 5.

**Murine Isotype:** IgG1, κ

**Specificity :** mouse and human keratin 5.

**Applications :**

- ELISA
- Western Blot
- Immunocytochemistry
- Immunohistochemistry

**Recommended dilutions:** 1/500-1/5000

Storage : store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

Does not contain any preservative

