

## **Anti-cell membranes (Membrane marker) mAb**

**Ref:** AO-01-MB-3H12

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 5Tγ-3H12

**Immunogen:** Natural immunisation, Unknown immunogen.

**Murine Isotype:** IgM, κ

**Specificity:** human and mouse cell membranes. Detects internal and plasma cell membranes.

**Applications :**

- Immunocytochemistry

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles.

***Does not contain any preservative***

