

Monoclonal Anti-MALTOSE BINDING PROTEIN Anti-MBP

Ref: AO-01-MBP-AS

Volume: 10, 25, 50 and 100 μ l

Clone: 17TF2-1H4

Format : Crude Ascites Fluid

Immunogen: MBP fusion Protein

Murine Isotype: IgG1 κ

Specificity: MBP and MBP fusion proteins

Applications :

- Immunoblotting
- ELISA
- Immunoprecipitation! not tested

Recommended dilutions: 1/1000-1/5000

Storage: Store at -20°C (up to 3 years). Avoid repeat freezing and thawing cycles.

Does not contain any preservative