

## **Monoclonal Anti-N-CoR**

**Ref:** AO-01-NCoR-1D10

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 1NC-1D10

**Immunogen:** recombinant human N-CoR protein, C terminal region (aa 1629-2440)

**Murine Isotype:** IgG1κ

**Specificity:** human N-CoR protein. Epitope between aa 1629 and 2440

**Applications :**

- Immunoblotting
- ELISA
- Immunofluorescence

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years).

Avoid repeat freezing and thawing cycles.

*Does not contain any preservative*

