

Monoclonal Anti-acetylated (K12) Histone H4

Ref: AO-01- H4-4C10

Volume: 10, 25, 50 and 100 µl

Format: Crude ascites fluid

Clone: 3HH4-4C10

Immunogen: synthetic peptide: MSGRG(acetyl-K)GG(acetyl-K)GLG(acetyl-K)GGA(acetyl K)RHRKVLR-(C) corresponding to aa 1 to 24 of human histone H4. This sequence is conserved in a wide range of species.

Murine Isotype: IgG1, κ

Specificity: Acetylated (K12) Histone H4.

Applications :

- ELISA
- Western blot
- Immunoprecipitation

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years) Avoid repeat freezing and thawing cycles.

Does not contain any preservative