

## Monoclonal Anti-acetylated (K5 or K12) Histone H4

**Ref:** AO-01-H4-3E7

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 3HH4-3E7

**Immunogen:** synthetic peptide: MSGRG(acetyl-K)GG(acetyl-K)GLG(acetyl-K)GGA(acetyl K)RHRKVLR-(C) corresponding to aa 1 to 24 of human histone H4. This sequence is conserved in a wide range of species

**Murine Isotype:** IgG1, κ

**Specificity:** Acetylated (K5 or K12) histone H4..

**Applications :**

- ELISA
- Western blot
- Immunoprecipitation

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years) Avoid repeat freezing and thawing cycles.

***Does not contain any preservative***