

Monoclonal anti-dimethyl (K14) Histone H3

Ref: AO-01- H3-2B10

Volume: 10, 25, 50 and 100 µl

Clone: 5HH3-2B10

Format : Crude Ascites Fluid

Immunogen : Ovalbumin conjugated synthetic peptide: TARKSTGG[dimethyl-K] APR-(C) corresponding to aa 6-17 of Histone H3. This sequence is conserved in a wide range of species..

Murine Isotype: IgG1, K

Specificity! dimethylated dimethylated (K14) Histone H3 (17 kDa).

Applications :

- ELISA
- Western blot
- Immunoprecipitation

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years). Avoid repeat freezing and thawing cycles.

Does not contain any preservative