

## **Monoclonal Anti-dimethyl / trimethyl K64 Histone H3**

**Ref:** AO-01- H3-1F2

**Volume:** 10, 25, 50 and 100 µl

**Clone:** 10HH3-1F2

**Immunogen :** synthetic peptide coupled to ovalbumin (STELLIR[me3K]LPFQRLVG-C) corresponding to aa 57-72 of human Histone H3. This sequence is conserved in a wide range of species.

**Murine Isotype:**

**Specificity:** dimethyl or trimethyl Histone H3 (17 kDa).

**Applications :**

- ELISA
- Western blot
- Immunoprecipitation

**Recommended dilutions:** 1/500-1/5000

**Storage:** Store at -20°C (up to 3 years). Avoid repeat freezing and thawing cycles.

***Does not contain any preservative***