

Monoclonal Anti-mCherry

Ref: AO-01- DS-2C4

Volume: 10, 25, 50 and 100 µl

Clone: 1DS-2C4

Immunogen: mCherry recombinant protein

Murine Isotype: IgG1, κ

Specificity: GFP, GFP-fusion proteins

Applications:

- ELISA
- Western blot
- Immunofluorescence

Recommended dilutions: 1/500-1/5000

Storage: Store at -20°C (up to 3 years). Avoid repeat freezing and thawing cycles.

Does not contain any preservative

For Research only

