

HRP conjugated Goat Anti-Rabbit IgG (H+L)

Product Information

Product Name	Cat. No.	Spec.
HRP conjugated Goat Anti-Rabbit IgG (H+L)	AO-03-GB23303	100 µL

Description

HRP conjugated Goat Anti-Rabbit IgG (H+L) is used in Western blot, ELISA and IHC applications. HRP Horseradish Peroxidase, which can catalyze ECL reagents to produce chemiluminescence in Western detection, TMB produces blue in ELISA, and DAB produce brown precipitation in IHC or Western blot detection.

Storage and Handling Conditions

Store at -20°C, valid for one year.

Recommended dilution ratio

ELISA 1:5000-1:10000

WB 1:10000-1:20000

IHC 1:500-1:5000

Notes

1. Please adjust the dilution ratio according to the experimental results. Minimize repeated freezing and thawing.
2. For your safety and health, please wear experimental clothes and disposable gloves.

For Research Use Only!

