

Mini Hand-Held Homogenizer

IG-L13k



IGENE LABSERVE

Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

FEATURES

- ❖ Light and easy to operate, broad scope of application.
- ❖ Use micro-DC motor for grinding trace tissue, efficient and effective.
- ❖ Grinding pestle is light and reusable, suitable for various sterile ways, convenient to be Dismantled.



SPECIFICATIONS

Model	IG-L13K Homogenizer
Speed range	13000rpm
Working noise	50dB
Battery non-loaded working time	2 hours
Battery	AA battery x 2
Voltage	With AA batteries
Dimension (WxDxH)	38x32x200 mm
Net weight	0.18 kgs

