


# IGENE LABSERVE

Innovative • Interactive • Intuitive

## IG- 1008 CENTRIFUGE



Address: 16/2 |Ground Floor|  
Ashok Nagar| New Delhi-110018|  
Office: - 011-79613060|  
Email: [info@igenels.com](mailto:info@igenels.com)|  
Web: [www.igenels.com](http://www.igenels.com)



# IGENE LABSERVE

Innovative • Interactive • Intuitive

## Standard Features:

---

- ❖ Ideal for microfiltration and quick spin-downs
- ❖ Easy-to-use with dual start/stop function
- ❖ 7200rpm fixed speed
- ❖ Compatible with 2mL × 8 and PCR8 × 2 rotors
- ❖ Quiet running, noise level ≤ 45 Db
- ❖ Rotor change without tools
- ❖ Accelerates and brakes in seconds

## Technical Specifications:

---

Max. Speed	7200rpm
Max. RCF	2680×g
Rotor Capacity	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC110-240V/50Hz/60Hz
Noise level	≤45dB
Dimensions [D × W × H]	160 × 170 × 122mm
Weight	0.5kg

Address: 16/2 |Ground Floor|  
Ashok Nagar| New Delhi-110018|  
Office: - 011-79613060|  
Email: [info@igenels.com](mailto:info@igenels.com)|  
Web: [www.igenels.com](http://www.igenels.com)

